

MEDROMA

www.medroma.co.za

ORTHOTICS

PROSTHETICS

COMPRESSION GARMENTS



Leader in Phlebology

Product Catalogue 2023/2024






OUR HISTORY



Our history spans three centuries and four generations: constant evolution and reference point in the production of medical compression stockings and garments, distribution of products in phlebology, lymphology, vascular surgery, vulnology, aesthetic medicine and plastic surgery.





PRODUCTION PROCESS

Our production cycle is vertical and integrated: from selection, use and processing of high quality yarns to computer-controlled knitting machines, ending with precise finishing and final inspections of the products.



Raw materials
warehouse:
Control and
acceptance of
raw materials



Winding
department: yarn
treatment



Coiling
department:
covering of the
elastic yarn



Manufacture
and finishing
Department



Final control



Packaging



Finished
products
warehouse and
logistics

QUALITY CERTIFICATIONS

Quality has always been one of our most important objectives, from product certifications to company internal quality system. Our medical stockings are registered and approved by Organizations in many countries.



Gloria Stockings are Class 1
Medical devices and meet all
requirements of the EU Directive
93/42/EWG



Gloria stockings are FDA
approved



Stockings are routinely tested by
the Hohenstein Research Institute
in Germany



Gloria Compression stockings meet
the RAL-GZG 387 requirements for
German Quality Seal.



Gloria internal quality control
system is UNI EN ISO 9001:2015
and ISO 13485 certified.



Gloria stockings are ANVISA
approved in Brasil



GLORIA COMPRESSION THERAPY


Compression therapy is used in the treatment of venous diseases and lymphedema. The aim of treatment is to maintain a high quality way of life despite an insufficient venous system. The type and extent of treatment depends on the severity of the problem and a doctor should always be consulted.

Gloria medical stockings are prescribed by physicians throughout the world, come in a wide selection of fabrics, styles, sizes, colours and compression ranges.


HOW DOES COMPRESSION WORK?

Compression therapy works by applying controlled pressure to the surface veins, keeping their diameter small and forcing blood back into the deep vein system which in turn pushes the blood uphill towards the heart.

Stronger compression is applied at the ankle and decreasing up the leg. Improved blood flow aids in the healing of varicose veins, venous leg ulcer and lymphatic problems.



HEALTHY VEINS



VARICOSE VEINS

CONTINUOUS EVOLUTION

Continuous evolution of scientific research, progress in technology and ever-growing demands by patients dictate that a modern industrial group must continually strive for improvements. Participation to international congresses and trade fairs, constant ties with doctors, research institutions, orthopaedical engineers and medical experts are a critical factor to our success. This is what drives Gloria and why its products are appreciated all over the world.

THE BENEFIT OF GLORIA COMPRESSION THERAPY

- Great for tired, achy legs
- Helps reduce and prevent swelling
- Assists with venous return of blood flow
- Recommended for venous insufficiency and lymphedema
- Maintaining a good quality way of life

PREVENTION & COMPRESSION THERAPY

MEDICAL COMPRESSION STOCKINGS 7

LYMPHOEDEMA & LIPEDEMA
FLAT KNIT GARMENTS 17

EASYVEN 18

VENOUS ULCER 20

COMPRESSION ARM SLEEVES..... 22

MEDICAL SUPPORT STOCKINGS AND SOCKS 25

CONTACTS 32

SIZING CHART 33







MEDROMA

EVERYDAY COMFORT AND ELEGANCE ALSO **FOR MEN**

Medical compression socks made of thin microfiber yarn that combined with a modern ribbed design, highlights the male appearance and makes the product soft and elegant without losing the therapeutic effects of compression.

Available in three stylish colors that make these classic men's socks extremely versatile and comfortable.



MICROFIBER RIBBED DESIGN CL. 1 - 2



COMPRESSION CLASSES:

GENTLEMAN 199

Compression class 1

GENTLEMAN 299

Compression class 2

INDICATIONS:

CLASS 1 (18-21 mmHg)

Mild symptoms of venous diseases such as heaviness and fatigue in the legs, mild varicosities without significant edema, post sclerotherapy, prophylaxis of thrombosis.

CLASS 2 (23-32 mmHg)

Moderate chronic venous insufficiency (CVI), thrombophlebitis (superficial thrombosis) with or without edema, to guarantee the therapy success (post sclerotherapy treatment, ambulatory selective varicose vein ablation, stripping) prophylaxis of thrombosis and embolisms.

AVAILABLE STYLES

MEDROMA

GENTLEMAN
MICROFIBER RIBBED DESIGN CL. 1 - 2

CLOSED TOE



AD

CL.1	199
CL.2	299
COLORS	<div>●</div> <div>●</div> <div>●</div>


DETAILS

COMPOSITION:
Polyamide - Elastane

CERTIFIED MEDICAL STOCKINGS



COLORS



MACRO



SIZES (CM)

S - M - L - XL

NORMAL - MAX


SHORT - LONG

CIRCUMFERENCE MEASUREMENTS (CM)

	S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max
d	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55
b	20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (CM)

	SHORT	LONG
ℓD	34-38	39-44





MEDROMA

SHEERNESS, ELEGANCE AND MEDICAL COMPRESSION

Stockings and pantyhose suitable for women who look for softness and elegant appearance, with medically effective compression.

The use of state-of-the-art technologies and highly elastic yarns guarantee the best possible comfort and fit.

MICROFIBER CL. 1 - 2



COMPRESSION CLASSES:

MICRO 181/191

Compression class 1

MICRO 281/291

Compression class 2

INDICATIONS:

CLASS 1 (18-21 mmHg)

Mild symptoms of venous diseases such as heaviness and fatigue in the legs, combined venous and arterial insufficiency, mild varicosities without significant edema, post sclerotherapy, small varicose dilations of cutaneous veins, initial varices during pregnancy.

CLASS 2 (23-32 mmHg)

Moderate symptoms of venous diseases such as aching heaviness and fatigue in the legs, chronic venous insufficiency, post sclerotherapy, prophylaxis of thrombosis, moderate to severe varicosities, post phlebotic syndrome, lymphedema, stasis dermatitis due to venous insufficiency, post fracture or post surgical edema, edema from burns, varices during pregnancy, pregnant patients with previous phlebitis.


AVAILABLE STYLES

MEDROMA

MICRO
MICROFIBER CL. 1 - 2

CLOSED TOE


OPEN TOE



AD

AGH*

AT



AD

AGH*

AGG**

AT

CL.1	191	191	191	181	181	181	181
CL.2	291	291	291	281	281	281	281
COLORS							

* AGH with stay-up silicon band

** AGG: adjustable waist belt, left or right style.

SIZES (CM)

XS - S - M - L - XL

NORMAL - MAX


SHORT - LONG

CIRCUMFERENCE MEASUREMENTS (CM)

	XS		S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max	normal	max
cT	60-80	70-90	70-90	80-100	80-100	90-110	90-110	100-120	100-120	110-130
cH	→ 85	→ 95	→ 95	→ 105	→ 105	→ 115	→ 115	→ 125	→ 125	→ 135
cG	42-46	47-54	47-51	52-60	52-57	58-67	58-63	64-74	64-70	71-83
cF	37-40	41-47	41-45	46-52	46-50	51-59	51-56	57-65	57-62	63-73
cD	28-30	31-36	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55
cB	18-19		20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (CM)

	SHORT	LONG
lD	34-38	39-44
lF	52-59	60-68
lG	62-71	72-82



DETAILS

STAY-UP BAND:

Lace stay-up band

On demand: Stay-up band plain with silicon drops


COMPOSITION:

Polyamide - Elastane


CERTIFIED MEDICAL STOCKINGS



COLORS



MACRO





MEDROMA

GRADUATED MEDICAL COMPRESSION WITH ELEGANCE

Stockings and pantyhose where compression is combined with a delightful appearance. Fine yarn and availability of open or closed toe provide optimum comfort and satisfaction for patients.

Gloriamed Soft are particularly suitable for young and dynamic women who need a therapeutical effect without compromising elegance and style.

MULTIFIBER SYNTHETIC CL. 1 - 2



COMPRESSION CLASSES:

SOFT 141/151

Compression class 1

SOFT 241/251

Compression class 2

INDICATIONS:

CLASS 1 (18-21 mmHg)

Mild symptoms of venous diseases such as heaviness and fatigue in the legs, combined venous and arterial insufficiency, mild varicosities without significant edema, post sclerotherapy, small varicose dilations of cutaneous veins, initial varices during pregnancy.

CLASS 2 (23-32 mmHg)

Moderate symptoms of venous diseases such as aching heaviness and fatigue in the legs, chronic venous insufficiency, post sclerotherapy, prophylaxis of thrombosis, moderate to severe varicosities, post phlebitic syndrome, lymphedema, stasis dermatitis due to venous insufficiency, post fracture or post surgical edema, edema from burns, varices during pregnancy, pregnant patients with previous phlebitis.























AVAILABLE STYLES

MEDROMA

SOFT
MULTIFIBER SYNTHETIC CL. 1 - 2

CLOSED TOE

OPEN TOE

	AD	AG *	AGG **	AT	ATM	ATU		AD	AG *	AGG **	AT	ATM	ATU
CL.1	151	151	151	151	151	151		141	141	141	141	141	141
CL.2	251	251	251	251	251	251		241	241	241	241	241	241
COLORS													
													

• : Available color only for made to measure.
•• : Available styles ONLY custom made.

* AG: available also with silicon band.
** AGG: adjustable waist belt, left on right style.

DETAILS


STAY-UP BAND:
plain with silicon drops

COMPOSITION:
Polyamide - Elastane


CERTIFIED MEDICAL STOCKINGS



COLORS



MACRO



CIRCUMFERENCE MEASUREMENTS (CM)

	XS		S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max	normal	max
cT	60-80	70-90	70-90	80-100	80-100	90-110	90-110	100-120	100-120	110-130
cH	→ 85	→ 95	→ 95	→ 105	→ 105	→ 115	→ 115	→ 125	→ 125	→ 135
cG	42-46	47-54	47-51	52-60	52-57	58-67	58-63	64-74	64-70	71-83
cF	37-40	41-47	41-45	46-52	46-50	51-59	51-56	57-65	57-62	63-73
cD	28-30	31-36	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55
cB	18-19		20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (CM)

	SHORT	LONG
lD	34-38	39-44
lF	52-59	60-68
lG	62-71	72-82

SIZES (CM)


XS - S - M - L - XL

NORMAL - MAX

SHORT - LONG

NEW

CUSTOM MADE PRODUCTS AVAILABLE



FINE COTTON ELEGANCE AND COMFORT

A complete range of stockings where superior quality cotton, excellent smoothness and comfort are combined with a delightful look and therapeutical effect.

Hypoallergenic and transpiratory, they are well suited for those who suffer from dermatological problems. They are also ideal for warm and humid environments and for people with allergies to synthetic fibers.

FINE COTTON CL. 1 - 2



COMPRESSION CLASSES:

COMFORT 142/152

Compression class 1

COMFORT 242/252

Compression class 2

INDICATIONS:

CLASS 1 (18-21 mmHg)

Mild symptoms of venous diseases such as heaviness and fatigue in the legs, combined venous and arterial insufficiency, mild varicosities without significant edema, post sclerotherapy, small varicose dilations of cutaneous veins, initial varices during pregnancy.

























CLASS 2 (23-32 mmHg)

Moderate symptoms of venous diseases such as aching heaviness and fatigue in the legs, chronic venous insufficiency, post sclerotherapy, prophylaxis of thrombosis, moderate to severe varicosities, post phlebotic syndrome, lymphedema, stasis dermatitis due to venous insufficiency, post fracture or post surgical edema, edema from burns, varices during pregnancy, pregnant patients with previous phlebitis.

AVAILABLE STYLES

MEDROMA

COMFORT
FINE COTTON CL. 1-2

CLOSED TOE						OPEN TOE						
	AD	AGH	AT	ATM	ATU	AD	AF *	AG *	AGG**	AT	ATM	ATU
CL.1	152	152	152	152	152	142	-	142	142	142	142	142
CL.2	252	252	252	252	252	242	242	242	242	242	242	242
COLORS												
												

• AT model also available with back extra wide panty part.

* AF style available only in class 2 normal or with stay-up silicon band.

AG style: class 1 only with stay-up silicon band; class 2 normal or stay-up silicon band

** AGG: adjustable waist belt, left or right style.

DETAILS

STAY-UP BAND:

Plain with silicon drops

COMPOSITION:


Class 1: 38% Cotton
33% Polyamide - 29% Elastane

Class 2: 37% Cotton
34% Polyamide - 29% Elastane


CERTIFIED MEDICAL STOCKINGS



COLORS



MACRO



CIRCUMFERENCE MEASUREMENTS (CM)


	XS		S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max	normal	max
cT	60-80	70-90	70-90	80-100	80-100	90-110	90-110	100-120	100-120	110-130
cH	→ 85	→ 95	→ 95	→ 105	→ 105	→ 115	→ 115	→ 125	→ 125	→ 135
cG	42-46	47-54	47-51	52-60	52-57	58-67	58-63	64-74	64-70	71-83
cF	37-40	41-47	41-45	46-52	46-50	51-59	51-56	57-65	57-62	63-73
cD	28-30	31-36	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55
cB	18-19		20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (CM)

	SHORT	LONG
lD	34-38	39-44
lF	52-59	60-68
lG	62-71	72-82

SIZES (CM)

XS - S - M - L - XL
NORMAL - MAX
SHORT - LONG





MEDROMA

UNIQUE STRETCH CHARACTERISTICS

The optimum elasticity allows a constant compression over long periods of time. The special porous mesh assures a normal transpiration of the skin. A complete range of styles and sizes, including Custom-made stockings, fit most patient needs.

On demand only: all models are also available with natural rubber.

LATEX FREE CL. 2 - 3



COMPRESSION CLASSES:

STRONG 261

Compression class 2

STRONG 361

Compression class 3

INDICATIONS:

CLASS 2 (23-32 mmHg)

Moderate symptoms of venous diseases such as aching heaviness and fatigue in the legs, chronic venous insufficiency, post sclerotherapy, prophylaxis of thrombosis, moderate to severe varicosities, post phlebitic syndrome, lymphedema, stasis dermatitis due to venous insufficiency, post fracture or post surgical edema, edema from burns, varices during pregnancy, pregnant patients with previous phlebitis.

CLASS 3: (34-46 mmHg)

More severe presentations of indications for class 2; correctable lymphedema, severe chronic venous insufficiency, post phlebitic syndrome, chronic post surgical venous insufficiencies, severe post thrombotic syndrome, venous ulcers.

AVAILABLE STYLES

MEDROMA

STRONG
LATEX FREE CL. 2 - 3

OPEN TOE



	AD	AF*	AG*	AGG**	AT	ATM	ATU
CL.2	261	261	261	261	261	261	261
CL.3	361	361	361	361	361	361	361
COLORS							

DETAILS

STAY-UP BAND:

plain with silicon drops

COMPOSITION:

Polyamide - Elastam
(Latex Free)

CERTIFIED MEDICAL STOCKINGS




COLORS

34

Beige

MACRO



SIZES (CM)

XS - S - M - L - XL

NORMAL - MAX

SHORT - LONG


CUSTOM MADE PRODUCTS AVAILABLE

CIRCUMFERENCE MEASUREMENTS (CM)

	XS		S		M		L		XL		XXL
	normal	max	normal	max	normal	max	normal	max	normal	max	normal
cT	60-80	70-90	70-90	80-100	80-100	90-110	90-110	100-120	100-120	110-130	110-130
cH	→ 85	→ 95	→ 95	→ 105	→ 105	→ 115	→ 115	→ 125	→ 125	→ 135	→ 135
cG	42-46	47-54	47-51	52-60	52-57	58-67	58-63	64-74	64-70	71-83	71-78
cF	37-40	41-47	41-45	46-52	46-50	51-59	51-56	57-65	57-62	63-73	63-69
cD	28-30	31-36	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55	48-52
cB	18-19		20-22		23-25		26-28		29-31		32-34

LENGTH MEASUREMENTS (CM)

	SHORT	LONG
lD	34-38	39-44
lF	52-59	60-68
lG	62-71	72-82



FLAT KNIT COMPRESSION STOCKINGS FOR LYMPHOEDEMA & LIPOEDEMA

The special texture of flat-knit garments makes them suitable for compression therapy in phlebology, lymphology and scar therapy. The unyielding compression ability of flat-knit, together with the massage effect produced, encourages lymphatic drainage and reduces edema expansion. The flat and elastic linked seams fit the shape of the body perfectly and make them suitable for patients with either deformed limbs, extra large proportions or medical complications such as elephantiasis.



ABSOLUTE CONTRAINDICATIONS:

Advanced peripheral occlusive arterial disease (PAD), decompensated heart failure, septic phlebitis, acute deep thrombosis, phlegmasia coerulea dolens, acute erysipelas, acute bacterial, viral or allergic inflammation, untreated or advanced skin pathologies, contact dermatitis, untreated open wounds, lower and upper limbs edema of unknown cause, untreated malignant tumors.

RELATIVE CONTRAINDICATIONS:

Exudative dermatosis, intolerance to elastic stockings materials, sensitivity disorders of lower and upper limbs, presence of edema in advanced peripheral neuropathies, presence of edema in diabetes mellitus, chronic polyarthritis, malignant lymphedema, compensated heart failure, skin infections, skin diseases (e.g. dermatosis bullous). In presence of open wounds in the treated area, it is necessary to cover the wound area with a dressing or bandage before wearing the elastic compression garments.



**For all product details, consult
the dedicated catalog**

COMPRESSION CLASSES

STANDARD GARMENTS:

Class 1 - 2

CUSTOM MADE GARMENTS:

Class 1 - 2 - 3

INDICATIONS

Lymphoedema, Elephantiasis, Lipedema, Primary and secondary varicosis, Chronic venous insufficiency, Edematous swellings of various genesis, including post traumatic and postoperative edema, Edema and varicose changes in pregnancy, Thrombophlebitis and support in phlebothrombosis, Post thrombotic syndrome, Post vein surgery, after Sclerotherapy, Angiodysplasias, Dermatosclerosis, Congestive states as a result of immobility, Edema in presence of paresis or partial paresis of the lower extremities, Thrombosis prophylaxis, Prevention of venous ulcers.

MEDICAL COMPRESSION STOCKINGS IN SYNTHETIC

SOFT, LIGHTWEIGHT AND
EASY TO WEAR

Easyven[®] stockings, in soft and elastic synthetic yarns, combine fit and comfort without losing the therapeutic efficacy. Easy to put on and remove provide comfort and well-being throughout the day.



COMPRESSION CLASSES:

EASYVEN 1: 15 - 20 mmHg **NEW**

EASYVEN 2: 20 - 30 mmHg

INDICATION:

EASYVEN 1: In presence of light venous disorders such as heaviness and fatigue in the legs, even during pregnancy. In non-serious post-traumatic situations, post-plaster.

EASYVEN 2: moderate chronic venous insufficiency (CVI) and after: healing of superficial ulcer, phlebitis syndrome, sclerotherapy treatment, ambulant selective varicose surgery, stripping. Varices during pregnancy, prophylaxis of thrombosis/embolism.

CONTRAINDICATIONS:

Infections, eczema, ulcers of lower limbs of arterial origins, severe neuropathies, uncompensated heart failure, severe leg edema, pulmonary edema, pre-existing gangrene, severe arterial circulation disorders.

AVAILABLE STYLES















MEDROMA

EASYVEN


EASYVEN 1: 15 - 20 mmHg
EASYVEN 2: 20 - 30 mmHg

CLOSED TOE

OPEN TOE

	AD	AGH	AT	AD	AGH	AGG*	AT
CL.1	EASY 1	EASY 1	EASY 1	EASY 1	EASY 1	EASY 1	EASY 1
CL.2	EASY 2	EASY 2	EASY 2	EASY 2	EASY 2	EASY 2	EASY 2
COLORS							
							

* AGG: adjustable waist belt, left or right style.

 BLACK COLOR ONLY ON DEMAND

DETAILS

STAY-UP BAND:

plain with silicon drops

COMPOSITION:

Polyamide
Elastane

COLORS

34

Beige

60

Black

SIZES (CM)

S - M - L - XL

NORMAL - MAX


SHORT - LONG

CIRCUMFERENCE MEASUREMENTS (CM)

	S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max
g	47 - 51	52 - 60	52 - 57	58 - 67	58 - 63	64 - 74	64 - 70	71 - 83
d	31 - 34	35 - 40	35 - 38	39 - 44	39 - 42	43 - 50	43 - 47	48 - 55
b	20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (CM)

	SHORT	LONG
lD	34 - 38	39 - 44
lG	62 - 71	72 - 82



THE REVOLUTION IN THE TREATMENT OF VENOUS ULCERS



UlcerKit® is a therapeutic compression stocking system indicated for leg ulcers, easy to use improves patients' compliance. Gloria Ulcerkit® can also be used as a normal medical stocking when a certain effort for application is required.

COMPRESSION CLASSES:

UlcerKit® is available in **Class 1** and **2**. As the system is composed by two separated parts, in case of need of a higher compression near the ankle (where usually the ulcer is located), there is the possibility to combine the lower part in CL. 2 with the upper part in CL. 1.

So the range is the following:

- **UlcerKit® 20/20:** Ankle CL. 1, Calf CL. 1
- **UlcerKit® 30/20:** Ankle CL. 2, Calf CL. 1
- **UlcerKit® 30/30:** Ankle CL. 2, Calf CL. 2



POINT B



POINT C

CIRCUMFERENCE MEASUREMENTS (CM)

	S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max
cC	30-33	34-37	34-37	38-41	38-41	42-45	42-45	46-49
cB	19-21		22-24		25-27		28-30	

LENGTH MEASUREMENTS (CM)

ONE LENGTH

INDICATIONS:

Class 1: Treatment of venous ulcer associated with other compression aids. Chronic venous insufficiency, post phlebitic, post ulcer, after sclerotherapy. Prevention of recidivism. Initial varices during pregnancy.

Class 2/1: Treatment of venous ulcer when a higher ankle compression is required. Stronger action at the ankle and moderate at the calf.

Class 2: Same as class 1 and 1/2, but with a highest compression on the leg, edema.

CONTRAINDICATIONS:

Do not wear UlcerKit® in presence of severe arterial venous insufficiency (frequent in diabetics). In less severe cases UlcerKit® can be used and if not suitable for the patient the stocking can be removed.

Patented Product

IT 1319230 – EU 1337211

MULTI-LAYER STOCKING SYSTEM



Gloria UlcerKit[®] PRO is a multi-layer stocking system to treat venous ulcers and is available in two different compressions and fabrics:

- **UlcerKit[®] PRO 35 Cotton** (35–40 mmHg)
- **UlcerKit[®] PRO 40 Cotton** (40–45 mmHg)
- **UlcerKit[®] PRO 35 Strong** (35–40 mmHg)
- **UlcerKit[®] PRO 40 Strong** (40–45 mmHg)

PACKAGE CONTENTS:

UlcerKit[®] PRO 35:

Calf-length CL.2 (1 pcs)
Calf-length 18 mmHg (1 pcs)
Thin understocking AD (2 pcs)
Foot slip

UlcerKit[®] PRO 40:

Calf-length CL.2 (1 pcs)
Calf-length 24 mmHg (1 pcs)
Thin understocking AD (2 pcs)
Foot slip



POINT B



POINT C



POINT D

INDICATIONS:

UlcerKit[®] PRO 35

(35–40 mmHg):

Skin ulcers of the leg without venous insufficiency (eg post-traumatic, vasculitic), skin ulcers of the leg with mild to moderate venous insufficiency.

UlcerKit[®] PRO 40

(40–45 mmHg):

Skin ulcers of the leg with severe venous insufficiency.

CONTRAINDICATIONS:


Gloria UlcerKit[®] PRO should not be used in the presence of significant arterial disease with ankle-brachial index (ABI) of less than 0.8. Diabetes mellitus with diabetic neuropathy.

CIRCUMFERENCE MEASUREMENTS (CM)

	S		M		L		XL	
	normal	max	normal	max	normal	max	normal	max
cD	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55
cC	32-35	36-41	36-39	40-45	40-43	44-51	44-48	49-56
cB	20-22		23-25		26-28		29-31	

LENGTH MEASUREMENTS (CM)

	SHORT	LONG
lD	34-38	39-44





MEDROMA

ELASTIC COMPRESSION ARM-SLEEVES

Elastic compression GloriaMed® arm-sleeves lightweight and comfortable thanks to the use of fine synthetic yarns, provide excellent comfort to the superior limb.

FEATURES:

- Available styles: with or without mitten
- Standard sizes and custom-made
- 2 Lengths (short or long)
- Soft and comfortable to wear
- Thin and flattened seams



CONTRAINDICATIONS:

- Uncontrolled lymphedema
- Inflammation
- Open wounds
- Infections
- Arterial circulatory problems of the upper limb



COMPRESSION CLASSES:

GLORIA MED® SLEEVE 101:

Compression classes 1

GLORIA MED® SLEEVE 201:

Compression classes 2

INDICATIONS:

Arm-sleeves are indicated in the treatment of primary and secondary lymphedema of the upper limbs.

Class 1:

prevention of post-surgery lymphedema after mastectomy with or without chemo and radiological treatments. Treatment of mild lymphedema. Prophylaxis of recurrent edema after decongestive therapy.

Class 2:

Prevention of edema after decongestion is achieved. Treatment of secondary lymphedema. Post-mastectomy fluid accumulation and edema following radiation therapy. Post-surgery and post-traumatic edema after basic treatment.

AVAILABLE STYLES

MEDROMA

SLEEVE
SLEEVE 101 - 201





WITH SHOULDER AND STRAP

WITH OR WITHOUT SILICON BAND


COLORS

34

Beige

	CG	AG	CG	AG
CL.1	101	101	101	101
CL.2	201	201	201	201
COLORS				

All models are ambidextrous



SIZES (CM)

S - M - L - XL

SHORT - LONG

CUSTOM MADE PRODUCTS AVAILABLE

CIRCUMFERENCE MEASUREMENTS (CM)

	S	M	L	XL
cG	25-30	29-34	33-38	37-42
cD	21-25	24-28	27-31	30-34
cC	15-17	17-19	19-21	21-23

LENGTH MEASUREMENTS (CM)

	SHORT	LONG
C-G	41-46	46-52

FLAT KNIT SLEEVES AND GLOVES

IDEAL FOR THE TREATMENT OF UPPER LIMB LYMPHEDEMA

The special texture of flat-knit garments makes them suitable for compression therapy in Phlebology, Lymphology and Scar Therapy.

The unyielding compression ability of the flat-knit fabric, together with the produced massage effect, encourages lymphatic drainage and reduce edema expansion. The flat and elastic seams perfectly fit the body shape.

CHARACTERISTICS:

- Two-ways soft flat-knit fabric
- Flat and elastic seams
- Sleeves and gloves available in standard sizes and made to measure to ensure a perfect comfort.
- Long-lasting garments



COMPRESSION CLASSES

STANDARD GARMENTS:

Class 1 - 2

CUSTOM MADE GARMENTS:

Class 1 - 2 - 3

INDICATIONS:

Primary and secondary lymphedema of the upper limb (stage 2 and 3), lipoeedema (stage 2 and 3), post surgical and post traumatic edema, upper extremity deep vein thrombosis (Paget von Schrötter syndrome), hypertrophic scars (following burns for example), rheumatoid arthritis.

ABSOLUTE CONTRAINDICATIONS:

Ischemia (eg. Severe arterial insufficiency), decompensated heart failure, untreated septic phlebitis, ischemic thrombosis (Phlegmasia coerulea dolens). Rheumatic fibromyalgia or generally any form of intolerance due to prolonged pressure on the skin surface. Skin infections, acute or exudative dermatitis. Untreated malignant tumors.


RELATIVE CONTRAINDICATIONS:

Exudative dermatosis, intolerance to elastic garment materials, sensitivity disorders of upper limbs, presence of edema in advanced peripheral neuropathies, presence of edema in diabetes mellitus, chronic polyarthritis, malignant lymphedema, compensated heart failure, skin infections, skin diseases (e.g. dermatosis bullous).

In presence of open wounds in the treated area, it is necessary to cover the wound area with a dressing or bandage before wearing the elastic compression garments.

ELEGANCE AND MEDICAL COMPRESSION IN FASHIONABLE COLOURS

Attractive, sheer and comfortable stockings and pantyhose with graduated compression. Suitable for young and dynamic women who need a medical effect with elegance and style.



INDICATIONS:

8-10 mmHg

Support pantyhose physician-recommended for prolonged standing, sitting and for young ladies during pregnancy.

12-15 mmHg

Mild compression for tired and heavy legs. Prevention and support in presence of mild varicosity during pregnancy.

18-20 mmHg

Designed specially to provide therapeutic benefits and compression in presence of heaviness and fatigue in the legs, initial varices during pregnancy.

20-30 mmHg

Mild to moderate symptoms of venous diseases such as heaviness and fatigue in the legs, post-sclerotherapy and after surgery, mild to moderate varicosities without significant edema, in presence of more severe varices during pregnancy.

AVAILABLE STYLES

MEDROMA


LADY GLORIA

8-10 mmHg, 12-15 mmHg
18-20 mmHg, 20-30 mmHg

LADY GLORIA 8

8-10 mmHg - **Sheer** - Size 1-5

CLOSE TOE



AT*



ATM

NEW HONEY

SAHARA

DEER


SIERRA

BEAVER

BLU MARINE

DARK GREY

BLACK



* AT model has sizes 4 and 5 with extra wide gusset (panty).


● Running out model / Running out colors

LADY GLORIA 12

12-15 mmHg - **Fish net** - Size 1-5

CLOSE TOE

OPEN TOE



AD



AG



AGH




AT*



ATM


Stay-up
silicon band



AD



AGH



AT

Stay-up
silicon band

WHITE

NEW HONEY

SAHARA

DEER


SIERRA

BEAVER

BLU MARINE

DARK GREY

BLACK



* AT model has sizes 3, 4 and 5 with extra wide gusset (panty).

● Running out color


AVAILABLE STYLES

MEDROMA

LADY GLORIA 18

18-20 mmHg - **Fish net** - Size 1-5


CLOSE TOE



AD AG AGH AT* ATM

Stay-up
silicon band



OPEN TOE



AD AGH AT

Stay-up
silicon band

WHITE
NEW HONEY
SAHARA
DEER
SIERRA
BEAVER
BLU MARINE
DARK GREY
BLACK




* AT model has all sizes with extra wide gusset (panty).

LADY GLORIA 24

20-30 mmHg - **Fish net** - Size 1-6


CLOSE TOE



AD AGH AT*

Stay-up
silicon band


OPEN TOE



AD AGH AT

Stay-up
silicon band

NEW HONEY
DEER
BLACK



* AT model has all sizes with extra wide gusset (panty).

SIZE TABLE FOR LADY GLORIA 8/12/18 - AGH, AT, ATM STYLES

HEIGHT cm	WEIGHT Kg															
	43	46	49	52	55	58	61	64	67	70	73	77	81	85	88	
150	1															
153	1															
156			2													
159			2													
162						3										
165						3										
168																
171										4						
173										4						
176												5				
179												5				
181																
183																

SIZE TABLE FOR LADY GLORIA 24 - AGH, AT, ATM STYLES

HEIGHT cm	WEIGHT Kg																	
	43	46	49	52	55	58	61	64	67	70	73	77	81	85	88	92	95	99
150	1																	
153	1																	
156																		
159			2															
162																		
165						3												
168																		
171																		
173										4								
176												5						
179														5				
181																6		
183																		

SIZE TABLE FOR: AD STYLE ONLY

SIZES	SHOE SIZE	
	EU	USA
1	35 - 36	5 - 6
2	37 - 38	6½ - 7½
3	38 ½ - 39	8 - 8½
4	39 ½ - 40	9 - 9½
5	41 - 42	9½ - 10
6 *	42 ½ - 43	10½ - 11

DETAILS

COMPOSITION:
Polyamide - Elastane

* Size 6 available only for Lady Gloria 24



MEDROMA

MAN'SUPPORT[®]

S U P P O R T S O C K S

A COMPLETE RANGE FOR DYNAMIC MEN

Preventive socks of high quality, available in several colours and 2 compression selections. Their elegance and comfort permit wearing them throughout the day whether at work or leisure. The cotton style provides increased transpiration and reduced allergies.

In particular, the **18 Silver Cotton** line in addition to cotton contains a special yarn called **Q-SKIN[®]**. It is a polyamide fibre in which silver ions are inserted directly during the spinning process meaning that it limits bacterial growth preventing body odor. Thanks to antibacterial, odor-control, thermoregulation properties of silver the socks takes care of the skin, providing its wearer with breathability, freshness, hygiene and comfort.



INDICATIONS:

15-18 mmHg

Support socks physician-recommended for prolonged standing, sitting or travelling. For tired and heavy legs.

18-20 mmHg

Designed specially to provide therapeutic benefits and compression in presence of heaviness and fatigue in the legs. Prophylaxis of venous disorders in patients at risk.


AVAILABLE STYLES

MEDROMA

MAN'SUPPORT
15-18 mmHg, 18-20 mmHg


MAN'SUPPORT 15

15-18 mmHg



MAN'SUPPORT 18

18-20 mmHg



SIZES	SHOE SIZE	
	EU	USA
1	36 - 38	5 - 6
2	39 - 40	6½ - 7
3	41 - 42	8 - 8½
4	43 - 44	9½ - 10
5	45 - 46	11 - 11½
6 *	47 - 48	13 - 14

* Available Size only for ManSupport 18 Silver Cotton.
ManSupport 18 Cotton only for blue and black colors.

DETAILS

COMPOSITION:

Man'Support: Polyamide - Elastane

Man'Support Cotton: Polyamide - Cotton - Elastane

Man'Support Silver Cotton: Polyamide (Q-Skin®) Cotton - Elastane



COMPRESSION SOCKS FOR MAN & WOMAN

Sitting or standing for long periods of time without moving your legs is a condition that causes a constriction of the veins and a slow return of venous blood to the heart. Swollen and painful legs are the first possible signs.

The graduated compression socks for men and women "Travel & Work" contribute actively to the health of your legs.



AVAILABLE STYLES:

CLOSED TOE



AD

BLACK



BLUE



BEIGE



SIZES	SHOE SIZE	
	EU	USA
1	37 - 38	5 - 6
2	39 - 40	6½ - 7
3	41 - 42	8 - 8½
4	43 - 44	9½ - 10
5	45 - 46	11 - 11½

INDICATIONS:

15-20 mmHg

Ideal for long trips by car, plane, train, and for those who in everyday life maintains a static position for many hours, or with their legs crossed.

COMPOSITION:

Polyamide - Elastane





CONTACTS

Catherine Cassel
GENERAL MANAGER

MEDROMA

www.medroma.co.za

Cell: 083 274 7820
Office: 012 - 991 2227
Email: catherine@medroma.co.za
Orders: orders@medroma.co.za

ORTHOTICS

PROSTHETICS

COMPRESSION GARMENTS

Medroma (Pty) Ltd | Registration number: 2023/937651/07
Office 109 | First Floor | Soetdoring Office Park | Doringkloof | Centurion | 0157
PO Box 25 | Southdowns | Irene | 0123

Jolanda Hayes
BUSINESS DEVELOPMENT MANAGER

MEDROMA

www.medroma.co.za

Cell: 081 047 0531
Office: 012 - 991 2227
Email: jolanda@medroma.co.za
Orders: orders@medroma.co.za

ORTHOTICS

PROSTHETICS

COMPRESSION GARMENTS

Medroma (Pty) Ltd | Registration number: 2023/937651/07
Office 109 | First Floor | Soetdoring Office Park | Doringkloof | Centurion | 0157
PO Box 25 | Southdowns | Irene | 0123




Gloria Med S.p.A.

Via Diaz 7 - 22017 Menaggio (CO) - ITALY
Tel. (+39) 0344 32123 - Fax (+39) 0344 32035
e-mail: gloriamed@gloriamed.it
www.gloriamed.com



Circumference Measurements												
cm.	XS*		S		M		L		XL		XXL**	
	Normal	Max	Normal	Max	Normal	Max	Normal.	Max	Normal.	Max	Normal	Max
cG	42-46	47-54	47-51	52-60	52-57	58-67	58-63	64-74	64-70	71-83	71-78	79-86
cF	37-40	41-47	41-45	46-52	46-50	51-59	51-56	57-65	57-62	63-73	63-69	70-76
cD	28-30	31-36	31-34	35-40	35-38	39-44	39-42	43-50	43-47	48-55	48-52	53-58
cb	18-19		20-22		23-25		26-28		29-31		32-34	

Length Measurements		
cm.	Short	Long
A-D	34-38	39-44
A-F	52-59	60-68
A-G	62-71	72-82



* XS not available for Gloriamed Gentleman

** XXL available only for Gloriamed Strong 261/361 & Comfort 242