



MEDROMA

www.medroma.co.za

ORTHOTICS

PROSTHETICS

COMPRESSION GARMENTS



LYMPHOLOGY

Flat Knit garments



BURNS & SCARS

Compression garments

CONTACTS

Catherine Cassel

GENERAL MANAGER

Cell: 083 274 7820
Office: 012 - 991 2227
Email: catherine@medroma.co.za
Orders: orders@medroma.co.za

MEDROMA

www.medroma.co.za

ORTHOTICS

PROSTHETICS

COMPRESSION GARMENTS

Medroma (Pty) Ltd | Registration number: 2023/937651/07
Office 109 | First Floor | Soetdoring Office Park | Doringkloof | Centurion | 0157
PO Box 25 | Southdowns | Irene | 0123

Jolanda Hayes

BUSINESS DEVELOPMENT MANAGER

Cell: 081 047 0531
Office: 012 - 991 2227
Email: jolanda@medroma.co.za
Orders: orders@medroma.co.za

MEDROMA

www.medroma.co.za

ORTHOTICS

PROSTHETICS

COMPRESSION GARMENTS

Medroma (Pty) Ltd | Registration number: 2023/937651/07
Office 109 | First Floor | Soetdoring Office Park | Doringkloof | Centurion | 0157
PO Box 25 | Southdowns | Irene | 0123

www.medroma.co.za



Gloria Med S.p.A.

Via Diaz 7 - 22017 Menaggio (CO) - ITALY
C.F. / P.IVA IT 00230380131

**FLAT KNIT COMPRESSION
GARMENTS****4**

Lymphology	6
About Flat knit fabric	7
Indications and Contraindications	8

FLAT KNIT STANDARD

Upper limbs	10
Lower Limbs	13

FLAT KNIT CUSTOM MADE

Upper limbs	15
Lower Limbs	22

**EPITHELIZATION
GARMENTS****26**

Upper limbs	30
Lower Limbs	31
Masks	32
Body	33

**BURNS & SCARS
TREATMENT****34**

Masks and face bands	40
Body	42
Upper limbs	44
Lower Limbs	48
Compression Stump Shrinker	51

**CUSTOM MADE ORDER
FORMS****52**

Perfect Profile - Flat Knit garments	53
Perfect Skin - Scars & Burns garments ...	57






FLAT KNIT GARMENTS

for lymphoedema & lipoedema

Perfect Profile standard	10
Perfect Profile custom made	15

“Lymphoedema is a long-term (chronic) condition that causes swelling in the body’s tissues.”



LYMPHOEDEMA is the accumulation of excessive amounts of protein-rich fluid resulting in swelling of one or more regions of the body.

It is, therefore, a chronic degenerative and progressive pathology that manifests itself in the body with **swelling due to the accumulation of lymph** in the superficial connective tissues. This is due to a mechanical failure of the lymphatic system to completely reabsorb the lymph itself. The condition usually affects the limb(s) although it may also involve the trunk, breast, head and neck or genital area.

Lymphoedema requires constant treatment and early intervention is recommended to minimise time and age related changes. The swelling may progress without treatment. The skin is prone to thickening and the development of fibrosis and other secondary changes.

Left untreated, this stagnant, protein-rich fluid causes tissue channels to increase in size and number, reducing the availability of oxygen. This interferes with wound healing and provides a rich culture medium for bacterial growth that can result in infections: cellulitis, lymphangitis, lymphadenitis and skin ulcers.

Lymphoedema cannot be cured but it can be improved and managed with appropriate intervention (Education skin care, exercises, Manual Lymphatic Drainage, Compression Therapy).

The use of **COMPRESSION GARMENTS** is one of the cornerstones of lymphoedema treatment. They apply medically appropriate pressure to the swollen region. This means that the pooling of fluid can be reduced, the limb size and shape preserved and the lymphatic circulation supported and improved.

LIPDEMA, not to be confused with lymphoedema and obesity, is a rare disease that causes **excess fat to accumulate in the lower part of the body**.

Lipedema most often involves the buttocks, thighs and calves. The upper arms can also be affected. The condition does not affect the hands or feet.

The exact cause of lipedema is unknown. But the condition runs in families and may be inherited. This metabolic syndromic form mainly affects women, and usually starts or gets worse at the time of puberty, pregnancy or menopause. According to some studies 11% of the female population is affected by this disease even if the lack of diagnosis could lead us to believe that these data may be higher. **Compression therapy** is also essential for this pathology, in order to ensure a reduction of edema. The use of Flat Knit garments allow protection and containment and reduces the incidence of pain typical of this disease.

WHAT IS FLAT KNIT FABRIC?



A flat knit compression garment is distinguished by its technical characteristics, it is made on flat knitting machine which guarantees the correct compression and strength, ensuring a customized product, based on the patient's anatomy. Standard products are also available.

Thanks to the short-stretch properties of the knitted fabric, flat knit compression combines several positive features: the low resting pressure makes donning and doffing easier, whereas the high working pressure ensures increased compression during movement.

The seam of flat knit compression garments is typical, which is carried out paying particular attention to avoid unpleasant rubbing. A further focus is provided for those seams that follow the fingers and the areas of greatest sensitivity, giving comfort to the product for a better patient's compliance.

There is a wide range of products available both standard, for patients who have a regular anatomical shape, and made to measure.

Perfect Profile flat-knit products offer solutions for every patient's unique shape, delivering optimum compression, containment levels, daily comfort, essential for a correct and lasting treatment of these pathologies.

PRIMARY LYMPHOEDEMA

STAGE 0



STAGE 1



STAGE 2



STAGE 3



The special texture of **Perfect Profile** flat-knit garments makes them suitable for compression therapy in Phlebology, Lymphology and Scar therapy.

CHARACTERISTICS

- Two-ways soft flat-knit fabric
- Flat and elastic seams
- All garments are available in standard sizes and made to measure to ensure a perfect comfort.
- **STANDARD GARMENTS** compression classes: CL. 1 - 2
- **CUSTOM MADE GARMENTS** compression classes: CL. 1 - 2 - 3
- Long-lasting garments

Available colors:

beige  black  brown 

Composition:

Class 1 85% Polyamide - 15% Elastane
Class 2 80% Polyamide - 20% Elastane
Class 3 70% Polyamide - 30% Elastane

INDICATIONS

Lymphoedema, Elephantiasis, Lipedema, Primary and secondary varicosis, Chronic venous insufficiency, Edematous swellings of various genesis, including post traumatic and postoperative edema, Edema and varicose changes in pregnancy, Thrombophlebitis and support in phlebothrombosis, Post thrombotic syndrome, Post vein surgery, after Sclerotherapy, Angiodysplasias, Dermatosclerosis, Congestive states as a result of immobility, Edema in presence of paresis or partial paresis of the lower extremities, Thrombosis prophylaxis, Prevention of venous ulcers.

Primary and Secondary lymphedema of the upper limb (stage 2 and 3), lipoedema (stage 2 and 3), post surgical and post traumatic edema, upper extremity deep vein thrombosis (Paget von Schrötter syndrome), hypertrophic scars (e.g. following burns), rheumatoid arthritis.

CONTRAINDICATIONS

ABSOLUTE CONTRAINDICATIONS

Advanced peripheral occlusive arterial disease (PAD), decompensated heart failure, septic phlebitis, acute deep thrombosis, phlegmasia coerulea dolens, acute erysipelas, acute bacterial, viral or allergic inflammation, untreated or advanced skin pathologies, contact dermatitis, untreated open wounds, lower and upper limbs edema of unknown cause, untreated malignant tumors.

RELATIVE CONTRAINDICATIONS

Exudative dermatosis, intolerance to elastic stockings materials, sensitivity disorders of lower and upper limbs, presence of edema in advanced peripheral neuropathies, presence of edema in diabetes mellitus, chronic polyarthritis, malignant lymphedema, compensated heart failure, skin infections, skin diseases (e.g. dermatosis bullous). In presence of open wounds in the treated area, it is necessary to cover the wound area with a dressing or bandage before wearing the elastic compression garments.

WASH AND CARE GUIDE FOR FLAT- KNIT FABRIC:

it is recommended to hand wash the garments, use lukewarm water (do not exceed 40°C), with a mild soap and avoiding the use of fabric softeners. Do not iron, Do not dry clean. Dry flat avoiding the sun and other direct heat sources (heating or dryer). Do not bleach, if possible do not use oils, creams or ointments.




**LS-500****GAUNTLET WITH THUMB**

CLASS 1 - 2

PACK: 1 piece

CIRCUMFERENCE MEASUREMENTS (CM)

SIZE	1	2	3	4	5	6	
cA	16-17	18-19	20-21	22-23	24-25	26-28	
cC	15-16	16-17	17-18	18-19	19-20	20-21	



Available Colors


**LS-501N****GLOVE WITH FINGERS TO THE NAIL**

CLASS 1 - 2

PACK: 1 piece

CIRCUMFERENCE MEASUREMENTS (CM)

SIZE	1/14	1/15	1/16	2/15	2/16	2/17	3/16	3/17	3/18	4/17	4/18	4/19	5/18	5/19	5/20	6/19	6/20	6/21	
cC	14	15	16	15	16	17	16	17	18	17	18	19	18	19	20	19	20	21	
cA	16-17			18-19			20-21			22-23			24-25			26-28			



Available Colors





LS-503

ARM-SLEEVE NO MITTEN WITH STAY-UP SILICON BAND

CLASS 1 - 2

PACK: 1 piece

SLEEVE STYLE: AMBIDEXTROUS

STAY-UP BAND WITH SILICON DOTS: 3 cm high

Available Colors



LS-5031

ARM-SLEEVE WITH GAUNTLET AND STAY-UP SILICON BAND


CLASS 1 - 2

PACK: 1 piece

SLEEVE STYLE: AMBIDEXTROUS

STAY-UP BAND WITH SILICON DOTS: 3 cm high

Available Colors



CIRCUMFERENCE MEASUREMENTS (CM)

SIZE	1	2	3	4	5	6
cG	23-24	25-27	28-30	31-33	34-37	38-40
cC	15-16	16-17	17-18	18-19	19-20	20-21
cA	16-17	18-19	20-21	22-23	24-25	26-28

LENGTH MEASUREMENTS (CM)

	SHORT	LONG
lG	40,5-43,5	43,5-46,5

**LS-5032L****ARM-SLEEVE WITH LONG FINGERS (OPEN TOE) AND STAY-UP SILICON BAND**

CLASS 1 - 2

PACK: 1 piece

SLEEVE STYLE: AMBIDEXTROUS

STAY-UP BAND WITH SILICON DOTS: 3 cm high

Available Colors

**LS-505****ARM-SLEEVE NO MITTEN WITH SHOULDER AND WIDE SUPPORTING STRAP**


CLASS 1 - 2

PACK: 1 piece

SLEEVE STYLE: AMBIDEXTROUS

(front and back closure for right or left use)

Available Colors

**CIRCUMFERENCE MEASUREMENTS (CM)**

SIZE	1	2	3	4	5	6
cG	23-24	25-27	28-30	31-33	34-37	38-40
cC	15-16	16-17	17-18	18-19	19-20	20-21
cA	16-17	18-19	20-21	22-23	24-25	26-28

LENGTH MEASUREMENTS (CM)

	SHORT	LONG
lG	40,5-43,5	43,5-46,5

LS-5071

FOOT GLOVE WITH HEEL FASTENING

CLASS 1 - 2

PACK: 1 piece

MODEL: AMBIDEXTROUS WITH OPENING FOR THE SMALL TOE


N.B.: Standard foot glove supplied with heel fastening sewn on one side only to allow adjustment according to the length of the foot (detached side to be sewn).

Available Colors



CIRCUMFERENCE MEASUREMENTS (CM)

	1	2	3	4	5	6	7
cA1	18,5-19,5	20-21	21,5-22,5	23-24	24,5-25,5	26-27,5	28-29,5



**LS-508*****BELOW-KNEE STOCKING WITH STAY-UP SILICON BAND 3 CM HIGH**

CLASS 1 - 2

PACK: 1 piece

TOE: OPEN

Available Colors

**LS-513*****THIGH-LENGTH STOCKING WITH STAY-UP SILICON BAND 5 CM HIGH**

CLASS 1 - 2

PACK: 1 piece

TOE: OPEN

Available Colors




***Both styles (standard) are only available with stay-up silicon band and open toe.**

CIRCUMFERENCE MEASUREMENTS (CM)

SIZE	cA	cY	cB	cC	cE	cG
1	20-22	27-29	18-20	29-31,5	32-34,5	48-51
2	22-24	29-31	20-22	31,5-34,5	34,5-37,5	51-55
3	24-26	31-33	22-24	34,5-37,5	37,5-40,5	55-60,5
4	26-28	33-35	24-26	37,5-40,5	40,5-43	60,5-64
5	28-30	35-37	26-28	40,5-43,5	43-46,5	64-69
6	30-32	37-39	28-30	43,5-46,5	46,5-49,5	69-73
7	32-34	39-41	30-32	46,5-49,5	49,5-53	73-76
8	34-36	41-43	32-34	49,5-52,5	53-56,5	76-79

LENGTH MEASUREMENTS (CM)

	SHORT	REGULAR	LONG
lG	64-67	68-71	72-75
lD	34-36	37-39	40-42
lA	15,5-16,5	16,6-17,5	17,6-18,5



VARIATIONS FOR GLOVES SEE PAGE 21



LC-500

GAUNTLET WITH THUMB

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Colors



LC-501

GLOVE WITH FINGERS TO THE NAIL

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Gloves with long fingers are only available with open toe

Colors



LC-502

GLOVE WITH SHORT FINGERS (UP TO THE KNUCKLE)

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Colors





LC-5001

LONG GAUNTLET WITH THUMB

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Colors



LC-5011

LONG GLOVE WITH FINGERS TO THE NAIL

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Gloves with long fingers are only available with open toe

Colors



LC-5021

GLOVE WITH SHORT FINGERS (UP TO THE KNUCKLE)

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Colors



VARIATIONS FOR ARM-SLEEVES SEE PAGE 21



LC-503

ARM-SLEEVE NO MITTEN WITH STAY-UP SILICON BAND

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

STAY-UP BAND WITH SILICON DOTS: 3 cm or 5 cm high

Also available without stay-up silicon band (normal finishing)

Colors



LC-5031

ARM-SLEEVE WITH GAUNTLET AND STAY-UP SILICON BAND

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

STAY-UP BAND WITH SILICON DOTS: 3 cm or 5 cm high

Also available without stay-up silicon band (normal finishing)

Colors



LC-5032

ARM-SLEEVE WITH SHORT FINGERS (UP TO THE KNUCKLE) AND STAY-UP SILICON BAND

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

STAY-UP BAND WITH SILICON DOTS: 3 cm or 5 cm high

Also available without stay-up silicon band (normal finishing)

Colors



LC-5032L

ARM-SLEEVE WITH LONG FINGERS (OPEN TOE) AND STAY-UP SILICON BAND

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

STAY-UP BAND WITH SILICON DOTS: 3 cm or 5 cm high

Also available without stay-up silicon band (normal finishing)

Only available with open toe

Colors





LC-504

ARM-SLEEVE WITH SHOULDER AND NARROW SUPPORTING STRAP

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Colors



LC-5041

ARM-SLEEVE WITH GAUNTLET, SHOULDER AND NARROW SUPPORTING STRAP

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Colors



LC-5042

ARM-SLEEVE WITH SHORT FINGERS (UP TO THE KNUCKLE), SHOULDER AND NARROW SUPPORTING STRAP

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Colors



LC-5042L

ARM-SLEEVE WITH LONG FINGERS (OPEN TOE), SHOULDER AND NARROW SUPPORTING STRAP

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Only available with open toe

Colors



VARIATIONS FOR ARM-SLEEVES SEE PAGE 21



LC-505

ARM-SLEEVE WITH SHOULDER AND WIDE SUPPORTING STRAP

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Colors



LC-5051

ARM-SLEEVE WITH GAUNTLET, SHOULDER AND WIDE SUPPORTING STRAP

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Colors



LC-5052

ARM-SLEEVE WITH SHORT FINGERS (UP TO THE KNUCKLE), SHOULDER AND WIDE SUPPORTING STRAP

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Colors



LC-5052L

ARM-SLEEVE WITH LONG FINGERS (OPEN TOE), SHOULDER AND WIDE SUPPORTING STRAP

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Only available with open toe

Colors





LC-506

ARM-SLEEVE WITH SHOULDER ATTACHED TO BRA

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Colors



LC-5061

ARM-SLEEVE WITH GAUNTLET, SHOULDER ATTACHED TO BRA

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Colors



LC-5062

ARM-SLEEVE WITH SHORT FINGERS (UP TO THE KNUCKLE), SHOULDER ATTACHED TO BRA

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Colors



LC-5062L

ARM-SLEEVE WITH LONG FINGERS (OPEN TOE), SHOULDER ATTACHED TO BRA

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Only available with open toe

Colors



VARIATIONS FOR MADE TO MEASURE ARM-SLEEVES

- Pad on the back of the hand
- Only pocket for pad
- Pocket + Pad on the back of the hand
- Lining fabric (inside)
- VERTICAL stopper in silicon (side or front)
- HORIZONTAL stopper in silicon (side or front)

VARIATIONS FOR MADE TO MEASURE GLOVES

- Pad on the back of the hand
- Only pocket for pad
- Pocket + Pad on the back of the hand



LC-507

FOOT GLOVE

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT

Foot glove only available with open toe

Colors



LC-5071

FOOT GLOVE WITH HEEL FASTENING

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT - RIGHT


Foot glove only available with open toe

Colors



VARIATIONS FOR FOOT GLOVE

- Hole for the small finger





LC-508

BELOW-KNEE STOCKING

CLASS 1 - 2 - 3

PACK: 1 piece

VARIATIONS:

- Normal finishing
- Elastic band 2,5 cm high
- Stay-up silicon band with dots, 3 cm or 5 cm high

Open or closed toe

Colors



LC-513

THIGH-LENGTH STOCKING

CLASS 1 - 2 - 3

PACK: 1 piece

VARIATIONS:

- Normal finishing
- Stay-up silicon band with dots, 5 cm high

Open or closed toe

Colors





LC-509

ONE-LEG STOCKING WITH ADJUSTABLE WAIST BELT (LEFT OR RIGHT)

CLASS 1 - 2 - 3

PACK: 1 piece

STYLE: LEFT or RIGHT

Open or closed toe

WAIST BELT VARIATIONS:

- For hips with a circumference equal to / greater than 100 cm. => waist belt 10 cm high
- For hips with a circumference of less than 100 cm. => waist belt 8 cm high
- Model for children => waist belt 6 cm high

Colors



LC-510E

ONE-LEG STOCKING WITH ADJUSTABLE WAIST BELT (LEFT + RIGHT)

CLASS 1 - 2 - 3

PACK: 1 pair

Open or closed toe

Colors



VARIATIONS FOR CUSTOM-MADE STOCKINGS

- VERTICAL stopper in silicon (side or front)
- HORIZONTAL stopper in silicon (side or front)
- Lining fabric (inside)
- Scrotal pouch

Retromalleolar Pads

- At the ankle(s) (inside - outside or both)
- Pocket for pad
- Pocket + foam pad (thickness 5 or 8 mm)
- Pocket + Polystyrene pad (thickness 2 cm.)
- Pocket + rubber pad (thickness 1 mm)



LC-510

TIGHTS

CLASS 1 - 2 - 3

PACK: 1 pair

Open or closed toe



PANTY PART VARIATIONS:

- Closed panty with adjustable waist belt (velcro)
- Closed panty with elastic waist belt 3 cm high
- Front open panty with hooks closure

CROTCH VARIATIONS:

- Closed crotch (normal)
- Compressive closed crotch
- Open crotch

ADDITIONAL MODELS OF PANTYHOSE:



LC-510S

ONE LEG & HALF
TIGHTS (L-R)

LC-510M

MATERNITY
TIGHTS

LC-510U

MAN TIGHTS

LC-510PED

CHILD TIGHTS

Colors







EPITHELIZATION GARMENTS

Perfect Skin 30

EPITHELIZATION GARMENTS

MEDROMA

Line of **standard products** designed for different parts of the body, in several sizes and can be used in the various stages of burn healing.

The garments can be worn directly on the skin or over a light bandage once the wounds have healed or after the first neo-epithelization.

The curved parts such as the top of the head, under the chin, the elbow and ankle are made with 3D Flex fabric to allow a close and better fit. The main fabric is very soft and comfortable.

Machine washable.


CUSTOMIZABLE

The fabric can be cut to size with no danger of fraying.


For example, there is the possibility to cut the material that covers the head, adjust the measurements or create openings for the ears.




Available Colors



blue



beige



blue and pistacho

Color combination only for children garments

EPITHELIZATION GARMENTS



VEST INT – 801

LEGGINGS INT – 802

GLOVE INT – 803

ARM SLEEVE INT – 804

BELOW KNEE STOCKING INT – 805

OPEN-FACED HELMET INT – 806



















MAXILLOFACIAL SUPPORT INT – 807



INT 804

TUBULAR ARM-SLEEVE
















Softer fabric in the elbow area to facilitate movement

Size	Arm circumference (top)	Color
1	16 - 17 cm	 
2	18 - 19 cm	 
3	20 - 21 cm	 
4	22 - 23 cm	 
5	24 - 25 cm	 
6	26 - 27 cm	 
7	28 - 29 cm	 
8	30 - 31 cm	 
9	32 - 33 cm	 
10	34 - 35 cm	 



INT 803


GLOVE OPEN TOE

Size	Palm circumference	Color
1	12 - 13 cm	 
2	14 - 15 cm	 
3	16 - 17 cm	 
4	18 - 19 cm	 
5	20 - 21 cm	 
6	22 - 23 cm	 
7	24 - 25 cm	 
8	26 - 27 cm	 



INT 802

LEGGINGS WITH OPEN CROTCH


Size	Thigh Circumference	Color
1	28 - 30 cm	
2	31 - 33 cm	
3	34 - 37 cm	
4	38 - 41 cm	
5	42 - 45 cm	
6	46 - 49 cm	 
7	50 - 53 cm	 
8	54 - 57 cm	 
9	58 - 61 cm	 
10	62 - 65 cm	 



INT 805

BELOW KNEE STOCKING







Softer fabric in the heel area to facilitate movement

Size	Ankle circumference	Color
1	21 - 23 cm	
2	24 - 27 cm	
3	28 - 30 cm	
4	31 - 33 cm	
5	34 - 36 cm	 
6	37 - 39 cm	 
7	40 - 42 cm	 
8	43 - 45 cm	 



INT 806







OPEN-FACED HELMET

Size		Color
1	Child	 
2	Teenager and woman	 
3	Man	 



INT 807





















MAXILLOFACIAL SUPPORT

Size		Color
1	Child	 
2	Teenager and woman	 
3	Man	 



INT 801

VEST FRONT OPENING WITH HOOKS

Size	Chest circumference	Color
1	45 - 50 cm	 
2	51 - 56 cm	 
3	57 - 62 cm	 
4	63 - 68 cm	 
5	69 - 74 cm	 
6	75 - 80 cm	 
7	81 - 86 cm	 
8	87 - 92 cm	 
9	93 - 98 cm	 
10	99 - 104 cm	 
11	105 - 110 cm	 
12	111 - 120 cm	 



ELASTIC COMPRESSION GARMENTS

*for burns and treatment
of hypertrophic scars and keloids*

Perfect Skin 40



Burned patients, when they leave the hospital, have an appreciable appearance, but after a few weeks the hypertrophic process of the scars starts, which has a pejorative evolution and can cause serious malformations. The soft and rosy scars that initially cover the wounds are not very evident and the patient begins to feel comforted, having come out of this traumatic experience, but subsequently the scars change, becoming rougher to the touch and changing their color, at first pink to a darker color. What appears as normal scar tissue, soft and elastic, will give rise to red, bulky and hard hypertrophic scars. The inexorable development of plaques and unsightly thickenings represent not only the evident deterioration of the tissue, which can lead to deformity and disability, but also create psychological stress, limiting the progress of the patient's social rehabilitation. Once the hypertrophic scars have completed their process, surgery is the only solution to restore natural appearance and functionality to the skin.

COMPRESSIVE THERAPY

To reduce the formation of hypertrophic scars, the most used and effective method is the pressure exerted on the healing burn. Compression is performed with **elastic compression garments** (masks, vests, gloves, arm-sleeve, stockings) worn directly in contact with the area to be treated. In order for the compression garment to be effective, it must perfectly cover the anatomy of the area to be treated and for this reason **custom-made** products is often the best solution

WHY ELASTIC COMPRESSION THERAPY IN PRESENCE OF SCARS

A hypertrophic scar is a thick raised scar. It's an abnormal response to wound healing in which extra connective tissue forms within the original wound area. The result a raised scar. Normally, a small wound to the top layer of your skin heals nicely. New skin forms as the wound heals. With deeper wounds your body responds by making collagen to repair the wound. After three or four days the burn is infiltrated by fibroblasts which are responsible for the formation of scar collagen fibers, which stop proliferating after about three weeks, creating a balance between the production and distribution of collagen itself. Collagen is thicker than the rest of the skin. This thicker, less flexible tissue becomes a scar. Most scars are flat. However, sometimes the body makes extra collagen that results in a raised scar. This type of raised scar can be either a hypertrophic

**FIRST-DEGREE BURN****THIRD-DEGREE BURN**

scar or a keloid. If a wound closes within three weeks (spontaneously or after transplantation), the wound will be less hypertrophic, as it will contain fewer fibroblasts. The compression applied above the scar tissue counteracts the ultrafiltration of the arterial side of the capillaries, producing a relative hypoxemia of the scar, thus inhibiting the activity of the fibroblasts. The balance between production and degeneration of collagen in the scar is therefore positively influenced by the elastic compression. Compression therapy is widely used as a first line conventional treatment of problem scarring which is commonly seen in the case of more severe burn injuries, especially those that require a skin graft, or for post-operative scar complications.

The success of this therapy depends on the time of application (24 hours a day) for the entire maturation phase of the scar which goes from 6 to 12 months.

It is very important to educate and motivate the patient to use compression garments and to create perfectly customized elastic compression garments.

Wounds that respond positively to compression develop the following characteristics:

- **the scar becomes lighter because hyperemia decreases;**
- **the amount of scar tissue decreases and becomes softer.**

However, the total area of the scar remains unaffected.

So it is clear that elastic compression therapy is one of the main available techniques for the treatment of hypertrophic scars and keloids.

INDICATIONS

Management and prophylaxis of problem scars, including hypertrophic and keloid scars caused by skin damage following burns, trauma or complications of postoperative surgical scars.

CONTRAINDICATIONS

Uncompensated organ failure (eg heart, liver or kidneys), Deep vein thrombosis or untreated phlebitis, Severe arteriopathy (ABI 0.49 or less). In case of hypersensitivity towards one or more components of the product itself.

FABRIC:

- Soft to the touch
- Several colors available, depending on type of garments

DRESSMAKING:

- Custom made production and customization of each garment
- The few seams are exterior and completely flat to avoid rubbing in sensitive areas.
- Natural anatomic shape thanks to the particular knitting process.
- Openings can be made to help placement

PRESSURE:

- Pressure levels controlled at 22 to 50 mmHg
- Pressure levels of the garments adapted to patient's needs.
- Varying levels of pressure possible between different parts of the same garment.

COMPOSITION:

60% Polyamide - 40% Elastane

VARIATIONS TO THE MODELS AND ACCESSORIES:

For available variations or accessories, refer to the current catalogue and price list.

WASHING INSTRUCTIONS AND CARE:

It is recommended to wash the garment with neutral soap or mild detergent and at a maximum temperature of 40 °C.

Wash separately. Do not use fabric softener or detergent with the addition of softener. Rinse thoroughly, do not wring. Air Dry, away from direct heating or direct sun. Do not iron. Do not tumble dry. Avoid oil or grease stains, which result in garment damage.



SUGGESTIONS FOR USE

Regular medical examinations are required to ensure the success of the therapy. It is important that compression garments are worn on a regular basis for 12 to 24 months or until scar maturity and at least 23 hours a day. This guarantees maximum benefit for optimal results in scar management, both aesthetically and functionally. The garments are typically removed 3 times a day for bathing, cleaning the scar, and masks can be removed during meals. Once the elastic compression garment has been positioned, make sure it adheres perfectly to the skin and is snug but not painful.

Do not carry out repairs yourself but contact your local shop. It is possible that knots may form in the fiber: this is due to contact and rubbing with other clothes and does not represent a defect in the quality of the products.

WARNING: rings and other accessories as well as long nails can damage the garments. Toenails and calluses can cause damage as well as defective padding inside the shoes. To wear the garments it is therefore recommended to use seamless gloves.

WARRANTY:

ADULT GARMENTS:

6 months warranty, as long as washing and instructions of use are correctly followed

CHILDREN'S GARMENTS:

3 months warranty, as long as washing and instructions of use are correctly followed

Any garment manipulated by unauthorized personnel is exempt from warranty.

Available Colors



beige

black



pink

blue



brown

green

Color combination for children's garments



pink and mauve blue and pistacho

Gloves color combinations for adults and children



pink and mauve blue and pistacho Green and olive green Black and grey

The colors may vary from the pictures on the catalogue.



JM 018

ADJUSTABLE CHIN STRAP WITH COLLAR


Colors



JM 019 A

OPEN FACE MASK

Colors



JM 019 C

CLOSED FACE MASK WITH NASAL OPENING

Optional: nose piece

Colors




JM 021

FULL ADJUSTABLE CHIN STRAP WITH COLLAR

Colors






JM 022

**CHILD'S COLLAR
(0-14 YEARS)**

Colors



JM 024

ADJUSTABLE HEAD BAND

Colors



JM 025

FACE BAND WITH NASAL OPENING

Optional: nose piece

Colors



JM 026

FULL FRONT FACE MASK WITH NASAL OPENING

Optional: nose piece

Colors



HEAD

MEDROMA



JM 026 A

FACE MASK OPEN FOREHEAD

Optional: nose piece

Colors



JM 026 B

FACE MASK - CHINLESS

Optional: nose piece

Colors



BODY



JM 209

STERNUM SHAPER

Colors






JM 002

ADULT'S VEST

Variation: Front or rear opening, with hooks

Colors




JM 002 S

WOMAN'S SHAPED VEST

Variation: Front or rear opening, with hooks

Colors



JM 001

CHILD'S VEST (0-14 YEARS)

Variation: Front or rear opening, with hooks
Bottom edge with silicon strips

For babies/children is also available a BODY model. Front or back opening.

Colors



JM 003

JM 004

CHILD'S SLEEVELESS VEST (0-14 YEARS)

ADULT'S SLEEVELESS VEST

Variation: Front or rear opening, with hooks

Colors






JM 005

STAY-UP ARM-SLEEVE NO MITTEN

Silicon band 3 cm high

Colors



JM 005 2

STAY-UP ARM-SLEEVE WITH GAUNTLET

Silicon band 3 cm high

Colors



JM 007

JM 008

CHILD (0-14 YEARS)

ADULT

ARM-SLEEVE WITH SHOULDER STRAP, NO MITTEN

Colors





JM 005 1

SINGLE FINGER, OPEN OR CLOSED TOE

Optional: hand finger with wrist attachment

Colors



JM 006

GAUNTLET

Colors



JM 006 1

GAUNTLET WITH WEB SPACE STRIPS

Colors



JM 006 2



ELBOW-LENGTH GAUNTLET

Stay-up silicon band.

Also available with web space strips
(item code JM 00621)

Colors



**JM 016 1****GLOVE WITH WEB SPACE STRIPS****Toe:** open or closed**Colors****JM 017****ELBOW-LENGTH GLOVE WITH LONG FINGERS****Toe:** open or closed**Colors****JM 017 1****ELBOW-LENGTH GLOVE WITH LONG FINGERS AND WEB SPACE STRIPS****Toe:** open or closed**Colors**

N.B: All elbow-length gloves, or with a length of more than 10 cm from the wrist, have a 3 cm high stay-up silicon band (silicon with dots).

Please specify at the time of order if you do not require stay-up silicon band.




JM 016

GLOVE

Toe: open or closed

Colors






JM 005 3

FOOTLESS STOCKING

Stay-up silicon band 5 cm high

Colors




JM 009

CHILD'S FULL TIGHTS (0-14 YEARS)

Toe: open or closed

Panty: closure with hooks

Colors



JM 009 1

JM 010 1

CHILD'S BERMUDA

ADULT'S BERMUDA

Panty: closure with hooks

Colors



JM 010


ADULT'S TIGHTS

Toe: open or closed

Panty: closure with hooks

Colors






JM 011

ONE-LEG STOCKING (LEFT OR RIGHT)

Toe: open or closed

Also available the combination R+L with one single waist belt (adjustable).

Colors




JM 012

BELOW-KNEE STOCKING

Toe: open or closed

Top band: normal finishing or stay-up silicon band 3 cm high

Colors



JM 014

THIGH-LENGTH STOCKING WITH STAY-UP SILICON BAND

Toe: open or closed

Stay-up silicon band 5 cm high

Colors





JM 013

ANKLE SOCK

Toe: open

Colors



JM 016 P

FIVE TOE ANKLE SOCK

Toe: closed

Colors





JM 304

STOCKING FOR STUMP, HIP STRAP



JM 304 1

STAY-UP STOCKING FOR STUMP (L/R)

Stay-up band with silicon dots

JM 304 3

STOCKING FOR STUMP +
ONE SHORT LEG TO KNEE

JM 304 2

STAY-UP STOCKING FOR STUMP (L+R)

Stay-up band with silicon dots

JM 305

ARM-SLEEVE STUMP WITH SHOULDER
CAP AND VELCRO STRAP (L/R)

JM 305 1

STAY-UP ARM-SLEEVE STUMP (L/R)

Stay-up band with silicon dots

Colors



Perfect Profile flat knit garments	53
Perfect Skin burns and scars garments	57

Date New order YES ☐ NO ☐

Patient identification code/name

Age Man ☐ Woman ☐


Client / Doctor

Client's stamp and Signature

By placing this order I accept that Custom made garments cannot be returned or replaced, and declare to accept the general conditions of supply included in the current price list.

Quantity	LEFT	RIGHT	COMPRESSION CLASSES	CLASS 1	CLASS 2	CLASS 3
----------	------	-------	---------------------	---------	---------	---------

HAND



Glove length

☐ SHORT Glove AC¹

☐ LONG Glove AE

Models

☐ Gauntlet with thumb

☐ Glove with fingers to the knuckle

☐ Glove with fingers to the nail (open toe)

Variations

☐ Pad on the back of the hand

☐ Only pocket for pad

☐ Pocket + Pad on the back of the hand

Colors

☐ Beige

☐ Brown

☐ Black

Fingers circumference:

	cX	cY	cZ
①	<input type="text"/>	<input type="text"/>	<input type="text"/>
②	<input type="text"/>	<input type="text"/>	<input type="text"/>
③	<input type="text"/>	<input type="text"/>	<input type="text"/>
④	<input type="text"/>	<input type="text"/>	<input type="text"/>
⑤	<input type="text"/>	<input type="text"/>	<input type="text"/>

Fingers length:

	lY	lZ
①	<input type="text"/>	<input type="text"/>
②	<input type="text"/>	<input type="text"/>
③	<input type="text"/>	<input type="text"/>
④	<input type="text"/>	<input type="text"/>
⑤	<input type="text"/>	<input type="text"/>

* For babies or children specify lC-C¹ length in the notes.

Notes / Special Requests

Date _____	New order YES <input type="checkbox"/> NO <input type="checkbox"/>	Client's stamp and Signature _____
Patient identification code/name _____		
Age _____	Man <input type="checkbox"/> Woman <input type="checkbox"/>	
Client / Doctor _____		

By placing this order I accept that Custom made garments cannot be returned or replaced, and declare to accept the general conditions of supply included in the current price list.

Quantity	LEFT <input type="text"/>	RIGHT <input type="text"/>	COMPRESSION CLASSES	CLASS 1 <input type="checkbox"/>	CLASS 2 <input type="checkbox"/>	CLASS 3 <input type="checkbox"/>
----------	---------------------------	----------------------------	---------------------	----------------------------------	----------------------------------	----------------------------------

Models

Arm sleeve to wrist



Arm sleeve with gauntlet*



Arm sleeve with short fingers*



Arm sleeve with long fingers* (open toe only)



* ☐ one-piece (arm and glove) ☐ Separate pieces (arm + glove)

N.B. Write hand measurements on the custom form for gloves.


Colors

☐ Beige ☐ Brown ☐ Black


Variations

- ☐ Normal finishing ☐ Elastic band 3 cm high
- ☐ Stay-up silicon band ☐ 3 cm ☐ 5 cm
- ☐ With shoulder attached to bra (with velcro)
- ☐ With shoulder and narrow supporting strap
- ☐ With shoulder and wide supporting strap
- ☐ Silicon stopper VERTICAL _____ x _____ cm ☐ Side ☐ Front
- ☐ Silicon stopper HORIZONTAL _____ x _____ cm ☐ Side ☐ Front
- ☐ Elbow protection Lining (inside)

Circumference Measurements:




Length Measurements:






Notes / Special Requests

Date _____	New order YES <input type="checkbox"/> NO <input type="checkbox"/>	Client's stamp and Signature _____
Patient identification code/name _____		
Age _____	Man <input type="checkbox"/> Woman <input type="checkbox"/>	
Client / Doctor _____		

By placing this order I accept that Custom made garments cannot be returned or replaced, and declare to accept the general conditions of supply included in the current price list.

LEFT LEG	CIRCUMFERENCE MEASUREMENTS	RIGHT LEG
cT		cT
cH		cH
cK		cK
cG		cG
cF		cF
cE		cE
cD		cD
cC		cC
cB¹		cB¹
cB		cB
cY		cY
cA¹		cA¹
cA		cA

Compression classes
<input type="checkbox"/> Left leg CL: <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3
<input type="checkbox"/> Right leg CL: <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3
Quantity
_____ piece _____ pair
Models
<input type="checkbox"/> AD below knee
<input type="checkbox"/> normal finishing <input type="checkbox"/> elastic band <input type="checkbox"/> stay-up silicon band 3 cm
<input type="checkbox"/> AG thigh length
<input type="checkbox"/> normal finishing <input type="checkbox"/> stay-up silicon band 5 cm
<input type="checkbox"/> AT Tights <input type="checkbox"/> ATM Maternity tights <input type="checkbox"/> ATU Man tights
<input type="checkbox"/> AT/E Half tights (long leg + short leg)
<input type="checkbox"/> Single leg with waist belt <input type="checkbox"/> right <input type="checkbox"/> left
<input type="checkbox"/> Single leg with waist belt R+L (one belt)
Colors:
<input type="checkbox"/> Beige <input type="checkbox"/> Brown <input type="checkbox"/> Black
Foot
<input type="checkbox"/> open toe <input type="checkbox"/> closed toe

LEFT LEG	LENGTH MEASUREMENTS	RIGHT LEG		
front _____ cm (waist - groin length) back _____ cm (waist - gluteal fold length)				
lT		lT		
lH		lH		
lK		lK		
lG		lG		
lF		lF		
lE		lE		
lD		lD		
lC		lC		
lB¹		lB¹		
lB		lB		
				

Variations
<input type="checkbox"/> Closed panty with adjustable waist belt (velcro)
<input type="checkbox"/> Closed panty with elastic waist belt 3 cm high
<input type="checkbox"/> Front open panty with hooks closure
Crotch variations:
<input type="checkbox"/> Closed crotch (normal) <input type="checkbox"/> Closed compressive crotch
<input type="checkbox"/> Open crotch <input type="checkbox"/> Scrotal pouch
Stopper
<input type="checkbox"/> Silicon stopper VERTICAL cm _____ x _____ <input type="checkbox"/> Side <input type="checkbox"/> Front
<input type="checkbox"/> Silicon stopper HORIZONTAL cm _____ x _____ <input type="checkbox"/> Side <input type="checkbox"/> Front
Retromalleolar pads
<input type="checkbox"/> Ankle pad <input type="checkbox"/> inside <input type="checkbox"/> outside <input type="checkbox"/> or both
<input type="checkbox"/> Protection lining at: _____ Length. cm _____ Width cm _____

Notes / Special Requests

All lengths, up to point "lG", must be taken internally and with the patient lying down.

Date _____ New order YES ☐ NO ☐

Client's stamp and Signature _____

Patient identification code/name _____

Age _____ Man ☐ Woman ☐

Client / Doctor _____

By placing this order I accept that Custom made garments cannot be returned or replaced, and declare to accept the general conditions of supply included in the current price list.



Quantity	LEFT	COMPRESSION CLASSES:	CLASS 1	CLASS 2	CLASS 3
	RIGHT		CLASS 1	CLASS 2	CLASS 3

LEFT FOOT

Measures	Toe 1	Toe 2	Toe 3	Toe 4	Toe 5
Circumference cZ					
Circumference cX					
Length ℓ					

RIGHT FOOT

Measures	Toe 1	Toe 2	Toe 3	Toe 4	Toe 5
Circumference cZ					
Circumference cX					
Length ℓ					



Models
<input type="checkbox"/> Foot glove
<input type="checkbox"/> Foot glove with heel fastening
Variations
<input type="checkbox"/> Hole for the small finger (toe 5) _____ cm

Colors
<input type="checkbox"/> Beige
<input type="checkbox"/> Brown
<input type="checkbox"/> Black

Models
<input type="checkbox"/> Foot glove
<input type="checkbox"/> Foot glove with heel fastening
Variations
<input type="checkbox"/> Hole for the small finger (toe 5) _____ cm

Notes / Special Requests

Date

New order YES ☐ NO ☐

Client's stamp and Signature


Patient identification code/name

Age

Man ☐ Woman ☐

Client / Doctor

By placing this order I accept that Custom made garments cannot be returned or replaced, and declare to accept the general conditions of supply included in the current price list.



Quantity

Models

- ☐ Child's collar
- ☐ Adult's collar
- ☐ Adjustable head band
- ☐ Adjustable chin strap with collar
- ☐ Full adjustable chin strap with collar
- ☐ Mask: open face
- ☐ Mask: closed face
 - ☐ with nose piece ☐ without nose piece
- ☐ Nose band
 - ☐ with nose piece ☐ without nose piece
- ☐ Full front mask
 - ☐ with nose piece ☐ without nose piece
- ☐ Face mask open forehead
 - ☐ with nose piece ☐ without nose piece
- ☐ Face mask chinless
 - ☐ with nose piece ☐ without nose piece

Variations and / or accessories

- ☐ With Nose bridge ☐ Without nose bridge
- ☐ Padding
- ☐ Patch of silicon gel cm 12,5 x 15

Colors:

- ☐ Beige ☐ Brown

IMPORTANT: MARK INJURIES ON THE DRAWING
Notes / Special Requests

Date _____ New order YES ☐ NO ☐

Client's stamp and Signature _____

Patient identification code/name _____

Age _____ Man ☐ Woman ☐

Client / Doctor _____

By placing this order I accept that Custom made garments cannot be returned or replaced, and declare to accept the general conditions of supply included in the current price list.

Quantity _____ right _____ left



Glove length

- ☐ Wrist-length glove
☐ Elbow-length glove (with stay-up silicon band)

Models

- ☐ Single finger toe: ☐ closed ☐ open
☐ Gauntlet
☐ Gauntlet with web space strips
☐ Glove long fingers (open toe)
☐ Glove long fingers (open toe) with web space strips
☐ Glove long fingers (closed toe)*
☐ Glove long fingers (closed toe)* with web space strips

* Fingertips:

- ☐ pointed ☐ square

Accessories

- ☐ Padding
☐ Patch of silicon gel cm 12,5 x 15

Colors:


for available colors see the catalog.

Notes / Special Requests


Date _____	New order YES <input type="checkbox"/> NO <input type="checkbox"/>	Client's stamp and Signature _____
Patient identification code/name _____		
Age _____	Man <input type="checkbox"/> Woman <input type="checkbox"/>	
Client / Doctor _____		

By placing this order I accept that Custom made garments cannot be returned or replaced, and declare to accept the general conditions of supply included in the current price list. _____

CIRCUMFERENCE AND LENGTH MEASUREMENTS



ADDITIONAL MEASURES FOR WOMAN'S SHAPED VEST



**IMPORTANT: MARK INJURIES
ON THE DRAWING**

Notes / Special Requests

Quantity

Vest: _____

Arm sleeve: ☐ right _____ ☐ left _____

PEDIATRIC Models (0-14 years)

OPENING		
<input type="checkbox"/> Child's vest	<input type="checkbox"/> front	<input type="checkbox"/> back
<input type="checkbox"/> Child's Body	<input type="checkbox"/> front	<input type="checkbox"/> back
<input type="checkbox"/> Child's sleeveless vest	<input type="checkbox"/> front	<input type="checkbox"/> back
<input type="checkbox"/> Child's sleeveless body	<input type="checkbox"/> front	<input type="checkbox"/> back
<input type="checkbox"/> Child's arm sleeve with shoulder strap		
<input type="checkbox"/> Arm stump garment		
<input type="checkbox"/> with shoulder strap		
<input type="checkbox"/> stay-up silicon band (3 cm high)		

ADULT Models

OPENING		
<input type="checkbox"/> Vest	<input type="checkbox"/> front	<input type="checkbox"/> back
<input type="checkbox"/> Woman's shaped vest	<input type="checkbox"/> front	<input type="checkbox"/> back
<input type="checkbox"/> Sleeveless vest	<input type="checkbox"/> front	<input type="checkbox"/> back
<input type="checkbox"/> Sternum shaper (chest measurement only)		
<input type="checkbox"/> Arm sleeve (wrist length)		
<input type="checkbox"/> with shoulder strap		
<input type="checkbox"/> stay-up silicon band (3 cm high)		
<input type="checkbox"/> Stay-up Arm sleeve with gauntlet		
<input type="checkbox"/> Arm stump garment		
<input type="checkbox"/> with shoulder strap		
<input type="checkbox"/> stay-up silicon band (3 cm high)		

Accessories

- ☐ Padding
- ☐ Patch of silicon gel cm 12,5 x 15

Colors:

for available colors see the catalog.

Date	New order YES <input type="checkbox"/> NO <input type="checkbox"/>	Client's stamp and Signature
Patient identification code/name		
Age	Man <input type="checkbox"/> Woman <input type="checkbox"/>	
Client / Doctor		


By placing this order I accept that Custom made garments cannot be returned or replaced, and declare to accept the general conditions of supply included in the current price list.

CIRCUMFERENCE MEASUREMENTS

WAIST cT _____ cm **HIPS** cH _____ cm

LEFT LEG

cG
cF
cE
cD
cC
cB¹
cB
cY
cA¹
cA



RIGHT LEG

cG
cF
cE
cD
cC
cB¹
cB
cY
cA¹
cA

IMPORTANT: MARK INJURIES ON THE DRAWING

Quantity

_____ right _____ left

Foot

☐ Open toe ☐ Closed toe

PEDIATRIC Models (0-14 years)

- ☐ Child's bermuda (panty closure with hooks)
- ☐ Child's tights (panty closure with hooks)
- ☐ Child's one-leg with hip support ☐ right ☐ left
- ☐ Child's one-leg with hip support R+L (one belt)

ADULT Models



- ☐ Below Knee ☐ normal finishing ☐ stay-up (3 cm band)
- ☐ Stay-up Thigh-length (silicon band 5 cm)
- ☐ Bermuda (panty closure with hooks)
- ☐ One-leg with hip support ☐ right ☐ left
- ☐ One-leg with hip support R+L (one belt)
- ☐ Tights (panty closure with hooks)

LENGTH MEASUREMENTS

back _____ cm (waist - gluteal fold length)


LEFT LEG

ℓG
ℓE
ℓD



RIGHT LEG

ℓG
ℓE
ℓD



Variations and / or accessories

- ☐ Closed crotch ☐ Open crotch
- ☐ Padding
- ☐ Patch of silicon gel cm 12,5 x 15

Colors:

for available colors see the catalog.

Notes / Special Requests


Date _____	New order YES <input type="checkbox"/> NO <input type="checkbox"/>	Client's stamp and Signature _____
Patient identification code/name _____		
Age _____	Man <input type="checkbox"/> Woman <input type="checkbox"/>	
Client / Doctor _____		

By placing this order I accept that Custom made garments cannot be returned or replaced, and declare to accept the general conditions of supply included in the current price list. _____


Quantity
<input type="checkbox"/> Left leg _____ piece (s) <input type="checkbox"/> Right leg _____ piece (s)

CIRCUMFERENCE MEASUREMENTS	
WAIST	cT _____
HIPS	cH _____

LEFT LEG
cG
cF
cE ¹
cE
cD
cC
cB ¹
cB
cA ¹
cA



RIGHT LEG
cG
cF
cE ¹
cE
cD
cC
cB ¹
cB
cA ¹
cA

LENGTH MEASUREMENTS		
<p>ℓZ: measure from the waist to the gluteal fold</p>  <p>ℓY: measure from the groin to the desired point.</p>		
ℓY _____ cm	ℓZ _____ cm	
Models		
<input type="checkbox"/> Stump stocking with stay-up silicon band		
<input type="checkbox"/> Stump stocking with waist belt		
Accessories		
<input type="checkbox"/> Padding		
<input type="checkbox"/> Patch of silicon gel cm 12,5 x 15		
Colors:		
<input type="checkbox"/> Beige	<input type="checkbox"/> Brown	<input type="checkbox"/> Black


IMPORTANT: MARK INJURIES ON THE DRAWING

Notes / Special Requests


Date _____	New order YES <input type="checkbox"/> NO <input type="checkbox"/>	Client's stamp and Signature _____
Patient identification code/name _____		
Age _____	Man <input type="checkbox"/> Woman <input type="checkbox"/>	
Client / Doctor _____		

By placing this order I accept that Custom made garments cannot be returned or replaced, and declare to accept the general conditions of supply included in the current price list.

Models	Quantity	Accessories	Colors:
<input type="checkbox"/> Ankle sock	Left _____ toe: <input type="checkbox"/> open <input type="checkbox"/> closed	<input type="checkbox"/> Padding	for available colors see the catalog.
<input type="checkbox"/> Five toes ankle sock	Right _____ toe: <input type="checkbox"/> open <input type="checkbox"/> closed	<input type="checkbox"/> Patch of silicon gel cm 12,5 x 15	



IMPORTANT: MARK INJURIES ON THE DRAWING



MEDROMA



Medroma (Pty) Ltd
Centurion, South Africa
012 - 991 2227

Gloria Med S.p.A.
Via Diaz 7 - 22017 Menaggio (CO) - ITALY